

## Product Datasheet

### VILL antibody (orb543043)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rabbit polyclonal antibody to VILL   |
| <b>Reactivity</b>          | Human, Mouse, Rat  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | ICC, IF, WB  |
| <b>Immunogen</b>           | Recombinant fusion protein of human VILL(NP_056957.3).   |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal,50% glycerol,pH7.3.  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks.<br>For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>MW</b>                  | 96kDa  |
| <b>Uniprot ID</b>          | <b>O15195</b>  |
| <b>Dilution Range</b>      | WB: 1:200-1:2000   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |