



Product Datasheet

NLRP1 antibody (orb543151)


Description	Rabbit polyclonal antibody to NLRP1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human NLRP1(NP_127497.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	160KDa
Uniprot ID	Q9C000
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of SH-SY5Y cells



Immunofluorescence analysis of PC-12 cells



Immunofluorescence analysis of K-562 cells