



Product Datasheet

SFRP5 antibody (orb543299)

Description	Rabbit polyclonal antibody to SFRP5
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human SFRP5(NP_003006.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35-38KDa
Uniprot ID	Q5T4F7
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse pancreas



Western blot analysis of extracts of Mouse pancreas