


Product Datasheet


ENO1 antibody (orb543372)

Description	Rabbit polyclonal antibody to ENO1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	Recombinant fusion protein of human ENO1(NP_001419.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47KDa
Uniprot ID	P06733
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of C6 cells

...



Immunofluorescence analysis of NIH-3T3

C...



Immunofluorescence analysis of C6 cells

...