



Product Datasheet

BMPR1B antibody (orb543657)


Description	Rabbit polyclonal antibody to BMPR1B
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Recombinant fusion protein of human BMPR1B (NP_001194.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	O00238
Entrez	658
Dilution Range	IHC-P: 1:50-1:100, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat ovar...



Immunohistochemical staining of human co...



Immunohistochemical staining of human ma...