


Product Datasheet

TRPM1 antibody (orb543769)

Description	Rabbit polyclonal antibody to TRPM1
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human TRPM1(NP_002411.3).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7Z4N2
Dilution Range	IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse re...