



Product Datasheet

BAZ2A antibody (orb543879)


Description	Rabbit polyclonal antibody to BAZ2A (TIP5)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human BAZ2A(NP_038477.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	280-310kDa
Uniprot ID	Q9UIF9
Dilution Range	IHC-P: 1:50-1:100, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS ce...



Immunohistochemical staining of mouse te...



Immunohistochemical staining of human CO...