

Product Datasheet

ERBB3 (Phospho-Y1328) antibody (orb544330)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to ERBB3 (Phospho-Y1328) |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ICC, IF, IHC, WB |
| Immunogen | A synthetic phosphorylated peptide around Y1328 of human ERBB3(NP_001973.2). |
| Preservatives | PBS with 0.05% proclin300,50% glycerol,pH7.3. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 185kDa |
| Uniprot ID | P21860 |
| Dilution Range | WB: 1:500-1:1000, IHC-P: 1:50-1:100, IF: 1:100-1:200 |
| Expiration Date | 12 months from date of receipt. |