

Product Datasheet

RAF1 (Phospho-S43) antibody (orb544339)

Description	Rabbit polyclonal antibody to RAF1 (Phospho-S43)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S43 of human Phospho-Raf1-S43 (NP_002871.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P04049
Entrez	5894
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.