

Product Datasheet

EIF2AK2 (Phospho-T451) antibody (orb544381)

Description	Rabbit polyclonal antibody to EIF2AK2 (Phospho-T451)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T451 of human PKR/EIF2AK2 (NP_001129123.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P19525
Entrez	5610
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.