

Product Datasheet

NFKBIA (Phospho-S32) antibody (orb544383)

Description Rabbit polyclonal antibody to NFKBIA (Phospho-S32)

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications IF, IHC, WB

Immunogen A synthetic phosphorylated peptide around S32 of human IκBα (NP_065390.1).

Preservatives Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG

Clonality Polyclonal

Purity Affinity purification


MW 39kDa

Uniprot ID [P25963](#)

Entrez [4792](#)

Dilution Range WB: 1:500-1:2000

Expiration Date 12 months from date of receipt.



Western blot analysis of extracts of HeLa...