

Product Datasheet

GRIN2B (Phospho-Y1474) antibody (orb544391)

Description	Rabbit polyclonal antibody to GRIN2B (Phospho-Y1474)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y1474 of human GRIN2B(NP_000825.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	190KDa
Uniprot ID	Q13224
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.