

Product Datasheet

PFKFB2 (Phospho-S483) antibody (orb544396)

Description	Rabbit polyclonal antibody to PFKFB2 (Phospho-S483)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S483 of human PFKFB2(NP_006203.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55KDa
Uniprot ID	O60825
Dilution Range	WB: 1:500-1:2000, IP: 1:50-1:100
Expiration Date	12 months from date of receipt.