

Product Datasheet

CBL (Phospho-Y774) antibody (orb544397)

Description	Rabbit polyclonal antibody to CBL (Phospho-Y774)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y774 of human CBL (NP_005179.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	110kDa
Uniprot ID	P22681
Entrez	867
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.