

Product Datasheet

CDK9 (Phospho-T186) antibody (orb544399)

Description	Rabbit polyclonal antibody to CDK9 (Phospho-T186)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T186 of human CDK9(NP_001252.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43KDa
Uniprot ID	P50750
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.