

Product Datasheet

RELA (Phospho-T435) antibody (orb544403)

Description	Rabbit polyclonal antibody to RELA (Phospho-T435)
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T435 of human RelA(NP_068810.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58kDa
Uniprot ID	Q04206
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.