

Product Datasheet

TP53 (Phospho-S392) antibody (orb544416)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to TP53 (Phospho-S392) |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | A synthetic phosphorylated peptide around S392 of human p53(NP_000537.3). |
| Preservatives | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 53kDa |
| Uniprot ID | P04637 |
| Dilution Range | WB: 1:500-1:2000 |
| Expiration Date | 12 months from date of receipt. |