

Product Datasheet

TSC2 (Phospho-S939) antibody (orb544420)

Description	Rabbit polyclonal antibody to TSC2 (Phospho-S939)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S939 of human TSC2 (NP_000539.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P49815
Entrez	7249
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.