

Product Datasheet

GJA1 (Phospho-S368) antibody (orb544427)

Description	Rabbit polyclonal antibody to GJA1 (Phospho-S368)
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	A synthetic phosphorylated peptide around S368 of human GJA1(NP_000156.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43kDa
Uniprot ID	P17302
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.