


Product Datasheet

PRKAA1 (Phospho-T183/PRKAA2-T172) antibody (orb544429)

Description	Rabbit polyclonal antibody to PRKAA1 (Phospho-T183/PRKAA2-T172)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding T183 of human PRKAA1
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	62kDa
Uniprot ID	Q13131
Entrez	55625563
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...