



Product Datasheet

HNF4 antibody (orb5453)


Description	HNF4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human HNF4 alpha (301-400/464 aa)
Target	HNF4A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	52 kDa
Uniprot ID	P41235
Dilution Range	WB=1:500-1000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of rat BRL-3A lysa...



Immunohistochemical analysis of formalin...



WB analysis of human Hcclm3 Cell lysates...