



Product Datasheet

CD109 antibody (orb545649)

Description	Mouse monoclonal antibody to CD109
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	WERI-RB-1 retinoblastoma cell line
Target	CD109
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q6YHK3
Entrez	135228
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human monocytes (red-fille...



Flow cytometry surface staining pattern ...