



Product Datasheet

CD143 antibody (orb545661)

Description	Mouse monoclonal antibody to CD143
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	dendritic cells
Target	CD143
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P12821
Entrez	1636
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of human stimulated monocytes...



Flow cytometry surface staining pattern ...