


Product Datasheet

MICA/MICB antibody (FITC) (orb545720)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody against MICA/MICB conjugated to FITC |
| Reactivity | Human |
| Conjugation | FITC |
| Tested Applications | FC |
| Immunogen | Transfected C1R cells expressing MICA |
| Target | MICA/MICB |
| Preservatives | Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide |
| Concentration | 1 mg/ml |
| Storage | Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze. |
| Note | For research use only |
| Application notes | Flow cytometry: Recommended dilution: 2 µg/ml. |
| Isotype | Mouse IgG2a |
| Clonality | Monoclonal |
| Purity | Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography. |
| Dilution Range | Flow cytometry: Recommended dilution: 2 µg/ml. |
| Expiration Date | 12 months from date of receipt. |



Separation of Jurkat cells stained using...