


Product Datasheet

Mouse Flt-3 Ligand Protein (orb545725)

Description	Mouse Flt-3 Ligand, His Tag is expressed from human 293 cells (HEK293). It contains AA Gly 27 - Arg ...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Flt-3 Ligand
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 20.2 kDa. The protein migrates as 25-33 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Mouse Flt-3 Ligand, His Tag (orb545725) is expressed from human 293 cells (HEK293). It contains AA Gly 27 - Arg 188 (Accession # P49772-1).
MW	20.2 kDa
Uniprot ID	P49772
NCBI	NP_038548.3
Expiration Date	6 months from date of receipt.



Immobilized
Human Flt-3,
Fc Tag (Cat.
No...