

Product Datasheet

Human GFR alpha-like Protein (orb545746)

Description	Human GFR alpha-like, His Tag is expressed from human 293 cells (HEK293). It contains AA Ser 19 - Gl...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	GFR alpha-like
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 39.7 kDa. The protein migrates as 45-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_997293.2
Purity	95%
Source	Human GFR alpha-like, His Tag (orb545746) is expressed from human 293 cells (HEK293). It contains AA Ser 19 - Glu 351 (Accession # Q6UXV0-1).
MW	39.7 kDa
Uniprot ID	Q6UXV0
NCBI	NP_997293.2
Expiration Date	6 months from date of receipt.