

## Product Datasheet

### Human B7-H6 / NCR3LG1 Protein (orb545764)

<b>Description</b>	Human B7-H6 Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Asp 25 - Ser...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	B7-H6
<b>Preservatives</b>	PBS, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-10×His
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 28.6 kDa. The protein migrates as 38-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	NP 001189368.1
<b>Purity</b>	90%
<b>Source</b>	Human B7-H6 Protein, His Tag (orb545764) is expressed from human 293 cells (HEK293). It contains AA Asp 25 - Ser 262 (Accession # Q68D85-1).
<b>MW</b>	28.6 kDa
<b>Uniprot ID</b>	<a href="#">Q68D85</a>
<b>NCBI</b>	<a href="#">001189368.1</a>
<b>Expiration Date</b>	6 months from date of receipt.