

Product Datasheet

Cynomolgus TIM-3 / HAVCR2 Protein (orb545768)

Description	Cynomolgus TIM-3, His Tag is expressed from human 293 cells (HEK293). It contains AA Ser 22 - Arg 20...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TIM-3
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 21.6 kDa. The protein migrates as 30-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	EHH54703.1
Purity	90%
Source	Cynomolgus TIM-3, His Tag (orb545768) is expressed from human 293 cells (HEK293). It contains AA Ser 22 - Arg 201 (Accession # EHH54703.1).
MW	21.6 kDa
NCBI	EHH54703.1
Expiration Date	6 months from date of receipt.