

Product Datasheet

Human FGFR3 protein (orb54587)

Description	Recombinant human FGFR3 protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-B2M-tagged
Note	For research use only
Application notes	N-terminal 6xHis-B2M-tagged: N-terminal 6xHis-B2M-tagged27-194AA: 397-806AAFull Length: Cytoplasmic Domain
Protein Sequence	RLRSPKKGLGSPTVHKISRFP LKRQVSLESNASMSSNTPLVRIARLSSGEGPTLA NVSELELPADPKWELSRARLT LGKPLGEGCFGQVMAEAIGIDKDRAAKPVTVA VKMLKDDATDKDLSDLVSEMEMMKMIGKHKNIINLLGACTQGGLYVLVEYAAK GNLREFLRARRPPGLDYSFDTCPPPEQLTFKDLVSCAYQVARGMEYLASQKCIH RDLAARNVLVTEDNVMKIADFLARDVHNL DYYKTTNGRLPVKWMAPALFD RVYTHQSDVWSFGVLLWEIFTLGGSPYPGIPVEELFKLLKEGHRMDK PANCTHD LYMIMRECWAAPSQRPTFKQLVEDLDRVLTVTSTDEYLDLSAPFEQYSPGGQD TPSSSSGDDSVFAHDLLPPAPPSSGGSR T
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	59.4 kDa
Uniprot ID	P22607
Expiration Date	6 months from date of receipt.