

Product Datasheet

Lactacystin (orb545880)

Description	Lactacystin
Target	Proteasome
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Formula	$C_{15}H_{24}N_2O_7$
Purity	> 98%
MW	376.37
CAS Number	133343-34-7
Solubility (25°C)	Soluble in DMSO (up to 20 mg/ml), in Water (up to 10 mg/ml), or in Ethanol (up to 1 mg/ml).
Expiration Date	12 months from date of receipt.