

Product Datasheet

NVP-BEZ235 (orb546001)

Description	NVP-BEZ235
Target	Kinase
Form/Appearance	powder
Storage	-20°C
Note	For research use only
Formula	$C_{30}H_{23}N_5O$
Purity	> 98%
MW	469.55
CAS Number	915019-65-7
Solubility (25°C)	Soluble in DMSO (up to 1 mg/ml with warming) or in DMF (up to 10 mg/ml, dissolves slowly, requires heating to 75°C).
Expiration Date	12 months from date of receipt.