

Product Datasheet

Tozasertib (orb546008)

| | |
|--------------------------|------------------------------------|
| Description | Tozasertib |
| Target | Kinase |
| Form/Appearance | powder |
| Storage | -20°C |
| Note | For research use only |
| Formula | $C_{23}H_{28}N_8O_5$ |
| Purity | > 98% |
| MW | 464.59 |
| CAS Number | 639089-54-6 |
| Solubility (25°C) | Soluble in DMSO (up to 100 mg/ml). |
| Expiration Date | 12 months from date of receipt. |