

Product Datasheet

Resiquimod (R848) (orb546107)

Description	Resiquimod (R848)
Target	TLR (Toll-like receptor)
Form/Appearance	powder
Storage	-20°C
Note	For research use only
Formula	$C_{17}H_{22}N_4O_2$
Purity	> 98%
MW	314-38
CAS Number	144875-48-9
Solubility (25°C)	Soluble in DMSO (up to 30 mg/ml) or in Ethanol (up to 15 mg/ml with warming)
Expiration Date	12 months from date of receipt.