


Product Datasheet

Pig porcine BMP-2 ELISA Kit (orb546421)

Description	Pig porcine BMP-2 ELISA Kit
Reactivity	Porcine
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: Q283-R396
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	bone tissue, cell culture supernatants and serum.
MW	44702 MW
Uniprot ID	C8C060
Sensitivity	2 pg/ml
Expiration Date	Please enquire.



Pig BMP-2
PicoKine
ELISA Kit
standard
CU...