


## Product Datasheet

### Bovine TGF Beta 1 ELISA Kit (orb546451)

<b>Description</b>	Bovine TGF Beta 1 ELISA Kit
<b>Reactivity</b>	Bovine
<b>Range</b>	15.6 pg/ml - 1,000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: CHO; Immunogen sequence: A279-S390
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Sample Types</b>	cell culture supernatants, serum, plasma (EDTA) and urine.
<b>MW</b>	44341 MW
<b>Uniprot ID</b>	<b>P18341</b>
<b>Sensitivity</b>	1 pg/ml
<b>Expiration Date</b>	Please enquire.



Bovine TGF beta 1 PicoKine ELISA Kit sta...