


Product Datasheet

Pig porcine TGF Beta 1 ELISA Kit (orb546453)

Description	Pig porcine TGF Beta 1 ELISA Kit
Reactivity	Porcine
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: A279-S390
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (EDTA) and urine.
MW	44341 MW
Uniprot ID	P07200
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Pig TGF beta 1 PicoKine ELISA Kit standa...