


Product Datasheet

Human FLT4/VEGFR3 ELISA Kit (orb546459)

Description	Human FLT4/VEGFR3 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Y25-I776
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P35916
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Human
FLT4/VEGFR3
PicoKine
ELISA Kit
Sta...