


Product Datasheet

Human IL-17F/Interleukin-17F ELISA Kit (orb546467)

Description	Human IL-17F/Interleukin-17F ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: R31-Q163
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin or EDTA)
Uniprot ID	Q96PD4
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human IL-17F
PicoKine
ELISA Kit
Standard...