


## Product Datasheet

### Pig porcine TGF-Beta 2 ELISA Kit (orb546485)

<b>Description</b>	Pig porcine TGF-Beta 2 ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	31.2 pg/ml - 2,000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: NS0; Immunogen sequence: A303-S414
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Sample Types</b>	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
<b>MW</b>	47748 MW
<b>Uniprot ID</b>	<b>P09858</b>
<b>Sensitivity</b>	10 pg/ml
<b>Expiration Date</b>	Please enquire.



Pig TGF-beta 2 PicoKine ELISA Kit standa...