


Product Datasheet

Mouse Myostatin/GDF8 ELISA Kit (orb546496)

Description	Mouse Myostatin/GDF8 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D268-S376
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:50, concentration in serum and plasma is around 20 ng/ml
Sample Types	cell culture supernatants, tissue, serum and plasma (heparin or EDTA)
Uniprot ID	O08689
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse Myostatin/GDF8 PicoKine ELISA Kit ...