


Product Datasheet

Human GAS6 ELISA Kit (orb546573)

| | |
|----------------------------|---|
| Description | Human GAS6 ELISA Kit |
| Reactivity | Human |
| Range | 62.5 pg/ml - 4,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: NS0; Immunogen sequence: A49-A678 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Application notes | Dilution ratio of 1:10, concentration in serum and plasma is around 8 ng/ml |
| Sample Types | cell culture supernatants, serum and plasma (heparin or EDTA) |
| Uniprot ID | Q14393 |
| Sensitivity | 10 pg/ml |
| Expiration Date | Please enquire. |



Human GAS6
PicoKine
ELISA Kit
Standard C...