


Product Datasheet

Human CD206/CLEC13D ELISA Kit (orb546629)

| | |
|----------------------------|---|
| Description | Human CD206/CLEC13D ELISA Kit |
| Reactivity | Human |
| Range | 0.78 ng/ml - 50 ng/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: NS0; Immunogen sequence: L19-K1383 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Application notes | Dilution ratio of 1:10, concentration in serum and plasma is around 100 ng/ml |
| Sample Types | cell culture supernatants, serum and plasma (heparin, EDTA, citrate). |
| Uniprot ID | P22897 |
| Sensitivity | 50 pg/ml |
| Expiration Date | Please enquire. |



Human CD206/CLEC13D PicoKine ELISA Kit S...