


Product Datasheet

Mouse Vitamin D BP/GC/Vitamin D Binding Protein ELISA Kit (orb546643)

Description	Mouse Vitamin D BP/GC/Vitamin D Binding Protein ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L17-S476
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:50000, concentration in serum and plasma is around 286 ug/ml
Sample Types	cell culture supernatants, lysates, serum, plasma (heparin) and urine.
Uniprot ID	P21614
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Mouse
Vitamin
DBP
PicoKine
ELISA Kit
Sta...