


Product Datasheet

Human PTH1R ELISA Kit (orb546667)

Description	Human PTH1R ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: Y23-M189
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	Q03431
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
PTH1R
PicoKine
ELISA Kit
Standard ...