


Product Datasheet

Human BMP-2 ELISA Kit (orb546890)

| | |
|----------------------------|---|
| Description | Human BMP-2 ELISA Kit |
| Reactivity | Human |
| Range | 31.2 pg/ml - 2,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: CHO; Immunogen sequence: Q283-R396 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Sample Types | bone tissue, cell culture supernatants and serum. |
| MW | 44702 MW |
| Uniprot ID | C8C060 |
| Sensitivity | 2 pg/ml |
| Expiration Date | Please enquire. |



Human BMP-2 PicoKine ELISA Kit standard ...