

## Product Datasheet

### Human CCL22/MDC ELISA Kit (orb546905)

<b>Description</b>	Human CCL22/MDC ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.2 pg/ml - 2,000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: E.coli; Immunogen sequence: G25-Q93
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Sample Types</b>	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
<b>MW</b>	10625 MW
<b>Uniprot ID</b>	<b>000626</b>
<b>Sensitivity</b>	1 pg/ml
<b>Expiration Date</b>	Please enquire.

