


Product Datasheet

Human PDGF-AB ELISA Kit (orb546907)

Description	Human PDGF-AB ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: (S87-T211) + (S82-T190)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	27283 MW
Uniprot ID	P04085, P01127
Sensitivity	2 pg/ml
Expiration Date	Please enquire.



Human PDGF-AB PicoKine ELISA Kit standar...