


Product Datasheet

Mouse Gp130/IL6ST ELISA Kit (orb546919)

Description	Mouse Gp130/IL6ST ELISA Kit
Reactivity	Mouse
Range	125 pg/ml - 8,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q23-E617
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	102451 MW
Uniprot ID	Q00560
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
Gp130/IL6ST
PicoKine
ELISA Kit
sta...