


Product Datasheet

Human Resistin ELISA Kit (orb546920)

Description	Human Resistin ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: K19-P108
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	11419 MW
Uniprot ID	Q9HD89
Sensitivity	3 pg/ml
Expiration Date	Please enquire.



Human Resistin PicoKine ELISA Kit standa...