


Product Datasheet

Human EPCAM/Trop1 ELISA Kit (orb546928)

| | |
|----------------------------|---|
| Description | Human EPCAM/Trop1 ELISA Kit |
| Reactivity | Human |
| Range | 15.6 pg/ml - 1,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: NS0; Immunogen sequence: Q24 - K265 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Sample Types | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA). |
| MW | 34932 MW |
| Uniprot ID | P16422 |
| Sensitivity | 10 pg/ml |
| Expiration Date | Please enquire. |



Human
EPCAM
PicoKine
ELISA Kit
standard ...